BT20

Self-Contained Biological Indicator.

For Steam sterilization processes.

Usage

Monitoring vacuum-assisted and gravity air-displacement Steam Sterilization processes at 121-135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017 and ISO 11138-3:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

Classification

Class 1, according to risk (ANMAT).

Characteristics

Polypropylene tube: 45.2 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Polypropylene cap: 13.0 mm high x 10.9 mm external diameter. Wall thickness: 0.7 mm.

Cap filter: medical grade paper, 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.5 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium: 0.5 – 0.7 ml, **purple** color.

Spore carrier: filter paper, 7.0 mm diameter.

10⁵ or 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Upon completion of the sterilization cycle, incubate the indicator for a maximum of 24 hours at temperature between 55-62 °C.

Positive growth results (failure in the sterilization process): culture medium turns to **yellow** after incubation between 55-62 °C for 24 hours (maximum).

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is informed at 132 °C and 135 °C.

Z-Value is informed.

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep in a dark place in its original box.

Transport conditions

Storage conditions should be strictly followed.

Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

2 years.

Packing

30 units or 100 units per box.

Packing information: product code and description, intended use, presentation, regulation, bacterial strain, storage conditions, manufacturer information and data on box's label.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label printed in black. 1.5 mm chemical indicator line, printed with Steam reactive ink (color change to brown). Product code, batch number, expiration date, intended use and name of organism.

On product's packing: product code, batch number, bacterial load, manufacture and expiration date, barcode and datamatrix code.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

Possible target markets

Healthcare and Industry.

Other important information

It is advisable to incubate at 60 °C in Bionova® IC10/20 incubator.

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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