mectron

medical technology

USE AND MAINTENANCE MANUAL

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EN

IT

DE




FR


ES

PIEZOSURGERY[®] touch

CE
0051

PIEZOSURGERY[®] touch

Copyright

© Mectron S.p.A. 2021. All rights reserved. No part of this document can be reproduced in any form without the written consent of the copyright owner.

SUMMARY

| | | |
|-----------|---|-----------|
| 1 | Introduction | 1 |
| 1.1 | Intended Use | 1 |
| 1.2 | Description of the Device | 2 |
| 1.2.1 | Patient Group Directions | 2 |
| 1.2.2 | Patient Selection Criteria | 2 |
| 1.2.3 | Indications for Use | 3 |
| 1.2.4 | Users | 3 |
| 1.3 | Disclaimer | 3 |
| 1.4 | Safety Precautions | 4 |
| 1.5 | Symbols | 6 |
| 2 | Identification Data | 8 |
| 2.1 | Device Identification Label | 8 |
| 2.2 | Handpiece Identification Data | 8 |
| 2.3 | Inserts Identification Data | 9 |
| 3 | Delivery | 10 |
| 3.1 | List of the Components of the PIEZOSURGERY® touch | 10 |
| 4 | Installation | 11 |
| 4.1 | Safety Requirements in the Installation Phase | 11 |
| 4.2 | Connecting the Accessories | 12 |
| 5 | Use | 16 |
| 5.1 | Switching the Device On and Off | 16 |
| 5.2 | Description of the Keyboard | 16 |
| 5.3 | Safety Requirements Before and During Use | 21 |
| 5.4 | Instructions for Use | 24 |
| 5.5 | Important Information on Inserts | 27 |
| 6 | Maintenance | 28 |
| 7 | Disposal Method and Precautions | 28 |
| 8 | Technical Specifications | 29 |
| 8.1 | Electromagnetic Compatibility IEC/EN 60601-1-2 | 31 |
| 8.1.1 | Guidance and Manufacturer's Declaration - Electromagnetic Emissions | 31 |
| 8.1.2 | Accessible Parts of the Casing | 32 |
| 8.1.3 | Guidance and Manufacturer's Declaration - Electromagnetic Immunity | 33 |
| 8.1.3.1 | Power Connection A.C. Input | 33 |
| 8.1.3.2 | Points of Contact with the Patient | 35 |
| 8.1.3.3 | Parts Accessible to the Input / Output Signals | 36 |
| 8.1.4 | Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device | 37 |
| 9 | Troubleshooting | 38 |
| 9.1 | Diagnostic System and Symbols on Keyboard | 38 |
| 9.2 | Troubleshooting Quick Guide | 40 |
| 9.3 | Fuses Replacement | 42 |
| 9.4 | Sending to an Authorized Mectron Service Centre | 43 |
| 10 | Warranty | 44 |

PAGE INTENTIONALLY LEFT BLANK

PAGE INTENTIONALLY LEFT BLANK

1 INTRODUCTION

Carefully read this manual before proceeding with installation operations, use, maintenance or other intervention operations on the device. Always keep this manual at hand.

Important: To avoid damages to persons or things, read with particular attention all the paragraphs "Safety Requirements" present in the manual.

According to the degree of severity the safety requirements are classified with the following indications:

⚠ WARNING: (Always referred to damage to persons)

⚠ CAUTION: (Referring to possible damages to things)

The aim of this manual is to make the operators aware of the safety regulations, installation

procedures, instructions for proper use and maintenance of the device and its accessories. Use of this manual for aims other than those strictly linked to the installation, use and maintenance of the device is prohibited.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is engaged in continuously updating its products with possible changes to the components of the apparatus.

In case you encounter discrepancies between the descriptions found in this manual and the equipment in your possession you can:

- check for any available updates in the *section MANUALS* of MECTRON *website*¹;
- ask clarifications to Your Dealer;
- contact MECTRON After Sales Service.

1.1 Intended Use

PIEZOSURGERY® touch is a piezoelectric ultrasonic device intended for oral bone surgery in the following applications:

- osteotomy and osteoplasty techniques;
- implantology;
- periodontal surgery;
- orthodontic surgery;
- endodontic surgery;
- maxillofacial surgery.

PIEZOSURGERY® touch, and its inserts, can operate as a piezoelectric scaler in the following dental applications:

- scaling: all the procedures for the removal of bacterial plaque deposits and supra-gingival, sub-gingival, and interdental calculi, and for the removal of stains;
- periodontology: periodontal treatment for the scaling and root-planing/ debridement without damage to the periodontium, including cleaning and irrigation of the periodontal pocket;

- endodontics: all treatments for root canals preparation, irrigation, filling, condensation of the gutta-percha and retrograde preparation;
- restoration and prosthesis: cavity preparation, removal of denture, condensation of the amalgam, finishing of the cervical margin, inlay/onlay preparation.

⚠ WARNING: The device must be used in a dental practice or ambulatory.

⚠ WARNING: Risk of explosion.

The device cannot function in places where there is an atmosphere saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

⚠ WARNING: Qualified and specialized personnel.

The device must be used exclusively by specialized personnel such as Surgeon. To use the device no special training is requested. The use of the device does not cause side effects if it is used correctly. An improper use might cause tissues heating.

¹ manuals.mectron.com

⚠ WARNING: Use the device only for the intended use. Failure to observe this precaution may cause serious injuries to the patient, the operator, and damages/breakdowns to the device.

⚠ WARNING: Carefully read and apply the recommendations given in this manual to avoid all risk to patient and/or user safety. Failure to do so may lead to serious injury to the patient and/or operator.

EN

1.2 Description of the Device

The user interface has been optimized with the PIEZOSURGERY® touch making all the functions readily available by integrating them in the touch keyboard.

The PIEZOSURGERY® touch is a device that uses ultrasonic piezoelectric technology to generate mechanical microvibrations of the inserts, to effectively cut mineralized tissues. This allows an efficient and safe cutting which preserves the integrity of the osteotomized surfaces.

The micrometric, ultrasonic vibrations of the inserts provide greater precision and a selective cutting action compared to

traditional methods such as drills or oscillating saws (which act with macrovibrations), therefore minimizing traumatic effect on soft tissues.

The cavitation effect of the irrigating solution helps to keep the operatory field blood-free. This provides an optimal intra-operative visual control thus increasing safety, even in areas that are anatomically most difficult to access.

This medical device can be used on any patient of any age, weight, height, gender and nationality.

1.2.1 Patient Group Directions

This medical device is designed to be used with the following patient population:

- Infants;
- Children;
- Adolescents;
- Adults;
- Elderly.

This medical device can be used on any patient of any age, weight, height, gender and nationality.

1.2.2 Patient Selection Criteria

The use of the device is not recommended in the following cases:

1. Patients with active implantable medical devices (for example: pacemakers, hearing aids and/or other electromagnetic prostheses) without the prior authorization of their doctor;
2. Women who are pregnant or breastfeeding, due to the restrictions associated with the possible use of medical solutions such as anaesthetics
3. Patients with allergies;
4. Patients with pathologies or clinical conditions for which it is not recommended to perform surgery or for which it may represent a contraindication according to the

treating physician. These conditions may include, but are not limited to: heart disease, diabetes, cirrhosis, HIV infection, pregnancy or breastfeeding, radiation therapy, chemotherapy, immuno-suppressive therapy, allergies and psychiatric disorders;

5. Patients with unsuitable treatment sites.

All models of "Piezosurgery®" bone surgery devices are intended for professional use only. Therefore, the user is the only person able to decide if and how to treat their patients.

1.2.3 Indications for Use

The use of the device is suitable for all the appropriate patients (see chapter above) for whom a bone surgery treatment is prescribed, by the treating physician, within the intended use of the device (see *Chapter 1.1 on page 1*).

1.2.4 Users

The device must be used only by specialised and properly trained personnel, such as a surgeon physician / dentist, able-bodied adults of any weight, age, height, gender and nationality.

1.3 Disclaimer

The manufacturer, Mectron, refuses all specific or implicit liability and shall not be held liable for any direct or indirect injury to people and/or damage to objects occurring following incorrect procedures connected with the use of the device and its accessories.

The manufacturer, Mectron, cannot be held specifically or implicitly liable for any type of injury to people and/or damages to objects caused by the user of the product and its accessories and which took place in the following cases:

1. Use or procedures different than those specified in the intended use of the product;
2. The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 8 on page 29;
3. The device is not used in compliance with all the instructions and precautions described in this manual;
4. The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions;
5. The assembly operations, extensions, adjustments, updates, and repairs on the device are performed by personnel not authorized by Mectron;
6. Improper use, mistreatments, and/or incorrect interventions;
7. Any and all attempts to tamper with or modify the device, under any circumstance;
8. Use of non-original Mectron inserts that entail a finite damage to the threading of the handpiece, thus compromising correct operation and causing risk of harm to the patient;
9. Use of non-original Mectron inserts, used in accordance to designed and tested settings of Mectron original inserts. The correct use of the settings is guaranteed only with original Mectron inserts;
10. Lack of stock materials (handpiece, inserts, wrenches) to be used in the event of device stop due to fault or of inconveniences.

1.4 Safety Precautions

⚠ WARNING: Risk of explosion.
The device cannot operate in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

ⓘ CAUTION: In the case that the end user, when operating in his or her own medical study or clinic, must subject the electro-medical equipment and systems to periodical inspections in order to adhere to imposed requirements, the test procedures that must be applied to electro-medical equipment and systems to evaluate safety must be performed in line with norm EN 62353 'Electro-medical devices - Periodical inspections and tests to be performed after repair interventions on electro-medical devices'. The interval for periodic checks, in the intended operating conditions and described in this "Use and Maintenance" manual, is one year.

⚠ WARNING: Checking device status before the treatment.

Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

ⓘ CAUTION: The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

ⓘ CAUTION: To avoid any risk of electric shock this device must be grounded.

⚠ WARNING: When the lid of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® touch must not be operated. The parts in movement may cause harm to the operator.

⚠ WARNING: All the accessories of the new or repaired devices are not sterile. All new and repaired accessories are supplied in non-sterile conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilization Manual.

⚠ WARNING: Infections control.
To ensure maximum safety of the patient and the operator, before using all the reusable parts and accessories, make sure you have first cleaned and sterilized them, following the instructions of the Cleaning and Sterilization Manual.

ⓘ CAUTION: Contraindications. Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ WARNING: Before each use, examine each part looking for any damages. If damages are found, do not use.

⚠ WARNING: Breakage and wear-out of the inserts. High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert.

Deformed or otherwise damaged inserts are susceptible to breakage during their use.

These inserts must never be used.

If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them.

The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent the patient from ingesting fragments of broken inserts.

When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out. Use of a worn-out insert reduces the cutting performance and can cause necrosis of the bone surface treated.

During the intervention, frequently check that the insert is intact, especially in its apical part. During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use. Do not exert excessive pressure on the inserts during their use.

⚠ WARNING: Only use original Mectron inserts, accessories, and spare parts.

⚠ WARNING: Contraindications. Do not use the PIEZOSURGERY® touch on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ CAUTION: Contraindications. Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ WARNING: Contraindications. Interference from other equipment. An electrical scalpel or other electro-surgical units near the PIEZOSURGERY® touch device may interfere with its correct operation.

⚠ WARNING: Contraindications. Interference with other equipment.

Although PIEZOSURGERY® touch complies with standard IEC 60601-1-2, it may interfere with other equipment located nearby.

PIEZOSURGERY® touch must not be used near, or stacked on other equipment. Install PIEZOSURGERY® touch at a safe distance from life-support equipment. If necessary, however, prior their use, please check and monitor that the device and all other equipment function correctly in that layout.

⚠ CAUTION: No modification of this equipment is allowed.

⚠ WARNING: Personal injury. Take care to ensure that cables do not hinder the free movement of personnel.

⚠ CAUTION: US Federal Law (of the United States of America) restricts this device to sale by or on the order of a physician.

⚠ WARNING: Do not activate the device if the handpiece is faulty, damaged or broken. Replace the handpiece immediately.

⚠ WARNING: If an adverse event and/or serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the manufacturer indicated on the product label.

⚠ WARNING: PIEZOSURGERY® touch is a device intended for bone surgery. However, prolonged contact and/or excessive force of the instrument tip on soft tissues should be avoided as this may cause thermal and/or blunt injury. Take particular care when sharp inserts are used. Prolonged mechanical action of a sharp insert may also result in the soft tissue being cut. In close proximity to soft tissues/nerves, it is recommended to complete the cut with a blunt tip diamond coated insert, to minimize the potential for soft tissue damage.

1.5 Symbols

| Symbol | Description | Symbol | Description |
|---|---|---|---|
|  | Device compliant with Regulation (EU) 2017/745. Notified body: IMQ S.p.A. |  Electrical Safety ES 60601-1 | Nemko Mark UL CSA compliance |
|  | Medical Device |  | Caution, read the Instructions For Use |
|  | Operating Instructions |  | Manufacturer |
|  | Date of manufacture |  | Serial Number |
|  | Lot Number |  | Product Number |
|  | Single-Use |  | Do not resterilize |
|  | Use-by date |  | Sterilized with Ethylene Oxide (EO) |
|  | Non-sterile |  | Can be sterilized in autoclave up to a maximum temperature of 135° C |
|  | Type "B" applied part as per EN 60601-1 recommendation |  | Equi-potentiality |
|  | Alternating Current |  | Connection of the foot pedal |
| I | Activation switch "on" | 0 | Activation switch "off" |
|  | Functional earth |  | Caution, risk of electric shock |
|  | Biohazard |  | The device and its accessories must not be disposed of or treated as solid urban wastes |

| Symbol | Description | Symbol | Description |
|---|---|---|---|
|  | General warning ^{a)} |  | Do not touch moving parts ^{b)} |
|  | Temperature limitation for transport and storage |  | Humidity limitation for transport and storage |
|  | Atmospheric pressure limitation for transport and storage |  | Do not use if package is damaged |
| IP22 | International Protection Code of the mechanical casing. | IP20 | International Protection Code of the mechanical casing. |
| QTY.1 | Quantity of items in the package: 1 | Rx Only | For US market, only. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed physician |

Table 1 – Symbols

a) The symbol is represented by a yellow warning triangle and a black graphical symbol.

b) The symbol is represented by a red circle-with-slash and a black graphical symbol.

2 IDENTIFICATION DATA


A precise description of the model and device serial number will facilitate the After-Sales Service to answer the inquiries quickly and effectively.

Always provide this information when contacting an Authorised Mectron Service Centre.

EN


2.1 Device Identification Label

Every device has an identification label that reports the main technical characteristics and the serial number. The identification label is located under the device. The complete technical specifications are reported in *Chapter 8 on page 29*.

A separate label shows further symbols and characteristics of the device. This identification label is placed under the device.

NOTE: The complete list of symbols and their description are shown in *Chapter 1.5 on page 6*.

2.2 Handpiece Identification Data

Each handpiece is laser-marked with (refer to *Figure 1 at page 9*):

- its serial number (SN);
- its product number (REF);
- its date of manufacturing;
- the device name (PIEZOSURGERY® white);
- the Mectron logo;
- HIBC 128 Data Matrix;
- CE, MD, Applied Part and Sterilizable symbols.

NOTE: The complete list of symbols and their description are shown in *Chapter 1.5 on page 6*.



Figure 1 – Handpiece Identification Data

2.3 Inserts Identification Data

The following data are laser-marked on each insert:

- the name of the insert (Ref. 1);
- the Mectron logo (Ref. 2);
- and the lot number to which the insert belongs (Ref. 3);

3 DELIVERY

3.1 List of the Components of the PIEZOSURGERY® touch

Refer to the front endpaper.

PIEZOSURGERY® touch consists of:

- A. Device core system;
- B. Peristaltic pump;
- C. Irrigation bag support rod;
- D. PIEZOSURGERY® touch handpiece fixed holder;
- E. PIEZOSURGERY® touch handpiece movable holder;
- F. Use and Maintenance Manual, Cleaning and Sterilization Manual;
- G. Power Cord;
- H. Foot Pedal with bracket, cable and connector;
- I. Touch keyboard protective film;
- J. PIEZOSURGERY® touch handpiece complete with cord and LED light front cone;

⚠ CAUTION: The handpiece and the cord cannot be separated.

- K. Metal front cone;
- L. Irrigation set (distributed by Mectron);
- M. Inserts/insert kit
- N. K8 Mectron torque wrench;

⚠ CAUTION: Some inserts require specific accessories and tools for being correctly used (for example kits and / or dedicated torque wrenches). Always read and follow the assembly and cleaning instructions supplied with the insert.

- O. Peristaltic pump tubing;
- P. Case.

Available accessories:

- Compatible inserts and insert kits (produced by Mectron);
- Bone grafting kit - Ref. O1 (distributed by Mectron);
- Compatible PIEZOSURGERY® handpieces (produced by Mectron);
- Compatible Mectron torque wrenches (produced by Mectron);
- Inserts holder (distributed by Mectron);
- Surgical tray - small (distributed by Mectron);

NOTE: Refer to the MECTRON website for the list of available and compatible accessories.

Available accessories for the thermosinfection:

- Adapter for inserts (produced by Mectron);
- Adapter for handpieces with cord (distributed by Mectron).

PIEZOSURGERY® touch consists of accessories that can be ordered separately. Refer to the "Packing List" included in your package, so that you may exactly know the quantity and type of the accessories supplied with the device you have purchased.

The packaging of the device cannot undergo strong impacts as contains electronic components, therefore the transport and the storage must be carried out with particular care. Do not stack more boxes in order to avoid damaging the underlying packaging. All material shipped by Mectron is controlled at the time of dispatch.

The device is shipped appropriately protected and packaged.

Upon receipt of the device, check for any possible damage caused during transport and in case any damage and/or defects is found,

complain to the transporter.
Keep the packaging in case there is a necessity to send any item to a MECTRON Authorized Service Centre and to store the device during long periods of inactivity.

⚠ WARNING: Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

4 INSTALLATION

The device must be installed in a suitable and convenient place for its use.

⚠ WARNING: The place where the device is installed must meet the requirements found in the *Chapter 4.1 on page 11*.

ready for use or may need to be enabled by typing an activation key.

In case Your device requires an activation key, the procedures to be followed may vary from country to country.

Always refer to your dealer for details.

PIEZOSURGERY® touch can be purchased

4.1 Safety Requirements in the Installation Phase

⚠ WARNING: Contraindications. Interference with other equipment.
Although PIEZOSURGERY® touch complies with standard IEC 60601-1-2, it may interfere with other equipment located nearby. PIEZOSURGERY® touch must not be used near, or stacked on other equipment. Install PIEZOSURGERY® touch at a safe distance from life-support equipment. If necessary, however, prior their use, please check and monitor that the device and all other equipment function correctly in that layout.

⚠ WARNING: Contraindications. Interference from other equipment. An electrical scalpel or other electro-surgical units near the PIEZOSURGERY® touch device may interfere with its correct operation.

⚠ CAUTION: The electrical system of the premises in which the device is installed and used must be compliant to the norms in force and to the relative electrical safety precautions.

⚠ CAUTION: To avoid any risk of electric shock this device must be grounded.

⚠ WARNING: Do not install the device in places where there is a risk of explosion.
The device cannot operate in environments where anesthetic or flammable mixtures are present.

⚠ WARNING: When the lid of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® touch must not be operated. The parts in movement may cause harm to the operator.

⚠ WARNING: Install the device in a place protected against collisions or against accidental sprays of water or liquids.

⚠ WARNING: Do not install the device above or near heat sources. Foresee adequate air circulation around the device when installing it. Leave adequate space, especially near the fan placed on the back part of the device.

⚠ CAUTION: Do not expose the device to direct sunlight or to sources of UV light.

⚠ CAUTION: The device can be transported, but it must be handled with care when it is displaced. Position the foot pedal on ground, so that it can only be activated intentionally by the operator.

⚠ **CAUTION:** Before connecting the handpiece cord to the device, make sure that the electrical contacts are perfectly dry. If need be, dry them with compressed air.

⚠ **CAUTION:** Position the device in way so that the power plug is always, easily reachable, since this plug is considered as a disconnecting means.

⚠ **WARNING:** Before using the Mectron PIEZOSURGERY® irrigation set, check the integrity of the sterile packaging and inspect the product to exclude the presence of any damage. Do not use the irrigation set if the packaging is open or damaged. The irrigation set loses sterility if the packaging is broken or damaged. If the packaging is damaged, proceed with the disposal of the set. Do not re-sterilize and reuse the set.


⚠ **CAUTION:** Do not allow the machine body or foot pedal to become wet. If liquid should enter the machine body or foot pedal, this may cause damage.

⚠ **CAUTION:** No modification of this equipment is allowed.

4.2 Connecting the Accessories


Insert the irrigation bag support rod in the dedicated hole;

1

Insert the handpiece fixed holder in its housing.

2


NOTE: The handpiece fixed holder can be placed in 4 different positions: right side, front right side, front left side and left side;

⚠ **WARNING:** The handpiece fixed holder must be used only to hold the Piezosurgery® handpiece.


Connect the peristaltic pump tubing to the connector located on the handpiece cord connection;

3

Insert the connector of the handpiece cord with the writing UP facing upward inside the connection in the front part of the PIEZOSURGERY® touch;


4

Fit the peristaltic pump tubing into the peristaltic pump proceeding as follows:


5

- Open completely the peristaltic pump cover;

- Insert the tube into the peristaltic pump;
- Close the peristaltic pump cover;

6


⚠ WARNING: When the cover of the peristaltic pump is open, the foot pedal of the PIEZOSURGERY® touch **MUST NOT** be activated. The parts in movement may cause harm to the operator.

Connect the peristaltic pump tubing end to the irrigation set;

⚠ WARNING: The irrigation set is provided in sterile packaging. Inspect the packaging to check its integrity. Whether it results damaged, DO NOT use the product and proceed with the proper disposal.


7

Insert the spike of the irrigation set into the irrigation bag used for the treatment (bag not supplied). Place the irrigation bag onto the appropriate stand;

⚠ WARNING: The irrigation bag holder stand must be used with bags of up to a maximum of 1000 ml capacity.

8



Connect the foot pedal to the back of the device in the socket marked with the symbol , using the foot pedal cable plug, ensuring it clicks into place;

NOTE: The foot pedal is provided with a bracket that, if positioned vertically, enables it to be moved to the most appropriate place for use, with no need to touch it with hands.


The bracket can also be positioned horizontally, if not in use;

9

Insert the power supply cord into its connector on the back of the device and then into the wall socket;

10


EN

Equipotential plug: The device is equipped with an additional equipotential plug located on the rear of the console.

This plug is in accordance with DIN 42801. Insert the connector of the equipotential cord (optional) to the plug on the rear of the device's console.


The purpose of additional potential equalization is to reduce differences of potential which can occur during operation between the device's body and conductive parts of other objects within the medical environment.

11

Clean the touch keyboard with a soft cloth with low fiber-release and position the sterile protective film on the surface of the clean touch keyboard.

12

⚠ WARNING: The sterile protective films are supplied inside sterile packages. Inspect the package to check its integrity. In case it is damaged, don't use it and dispose it properly.

5 USE

5.1 Switching the Device On and Off

Switching the Device On


To switch the device on, set the switch button located on the left side of the device to position "I", taking care not to press the foot pedal down during this stage.

4 symbols light up on the device (ref. U inside the cover) and progressively switch off. When this sequence is completed, the device is set to default parameters and it is ready to use.

Switching the Device Off

To switch the device off, set the switch button located on the left side of the device to position "0", taking care not to press the foot pedal down during this stage. The device turns off.


NOTE: every time the device is switched on, the default settings are applied: "function" ENDO, "irrigation" 3, "light" OFF.

5.2 Description of the Keyboard

TOUCH KEYBOARD


The user can configure the touch keyboard by simply touching on the touch keyboard. Depending on the selected setting, the electronic feedback system automatically adjusts the correct operating frequency.

FUNCTIONS (Ref. R front endpaper)

Depending on the type of surgery, it is possible to choose one of the 6 options available from the “function” list:

- **ENDO:** dedicated to endodontic surgery and to the detachment of the sinus membrane;
- **PERIO:** dedicated to periodontal surgery;
- **SPECIAL:** dedicated only to the inserts for osteotomies with thickness of 0.35 mm and for prosthetic applications;
- **CANCELLOUS:** dedicated to the cutting and removal of poorly mineralized bone;
- **CORTICAL:** dedicated to the cutting and removal of highly mineralized bone;
- **IMPLANT:** dedicated to bone perforation in the technique of the implant site preparation.

EN


FILLING AND CLEANING THE IRRIGATION CIRCUIT (Ref. T front endpaper)

The device is equipped with the “pump/clean” key which, depending on the mode of use, allows to perform the PUMP function or the CLEAN function.

The PUMP function can be used at the beginning of the treatment, to flush the entire irrigation line up to the insert, so that the surgery can be started with the necessary irrigation (see Chapter 5.4 on page 24).

The CLEAN function allows to perform a cleaning cycle of the irrigation circuit.

This function must be performed when you have finished using the device and before cleaning, disinfecting, and sterilizing all the parts (see Cleaning and Sterilization Manual).

IRRIGATION (Ref. Q front endpaper)

The capacity of the peristaltic pump can be adjusted by selecting the numbers on the touch keyboard in the "irrigation" column.



7 capacity levels are available:

- 0 = no operation of the pump: no irrigation outflows from the insert
- Da 1 a 6= the pump flow ranges from 8 ml/min to approximately 75 ml/min.

The irrigation levels available depends on the type of function selected, in the particular:


- ENDO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min)
- PERIO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min)
- SPECIAL - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)
- CANCELLOUS - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)
- CORTICAL - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)
- IMPLANT - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min)

NOTE: treatment without irrigation is possible only with the ENDO and PERIO functions, setting the irrigation capacity level on "0".

In the event that an irrigation capacity lower than 8 ml/min is required, use the "bone grafting kit" (accessory that can be ordered separately), inserting it between the irrigation set and the tube in silicon of the handpiece, making it pass through the peristaltic pump, and selecting 1 as the level of irrigation.

⚠ WARNING: If the "bone grafting kit" tube is kept inserted for the entire duration of the intervention, the capacity of the pump on all its levels is limited, independently of which insert is used.



LIGHT (Ref. S front endpaper)

Depending on the type of handpiece and operation that needs to be performed, it is possible to choose 3 options from the "light" list:

- By choosing the AUTO option, the LED light on the front terminal of the handpiece is switched on by pressing the foot pedal, and automatically switches off 3 seconds after the pedal is released.
- By choosing the ON option, the LED light on the front terminal of the handpiece stays on permanently, regardless of whether pressure is applied on the foot pedal or not. The light switches off 100 seconds after the last pressure of the foot pedal and the option shifts from ON to AUTO.
- By choosing the OFF option, the LED light on the front terminal of the handpiece stays off permanently.

The position of the LED light on the front terminal of the handpiece can be adjusted in the following way:


1. Hold the body of the handpiece and lightly unscrew the metal ring nut located at the base of the front terminal, rotating it counter-clockwise.
2. Rotate the front terminal so that the LED light goes into the desired and necessary position.
3. To fasten it into position, screw the metal ring nut, rotating it clockwise.

SYMBOLS (Ref. U front endpaper)

PIEZOSURGERY® touch is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol.

To help the user identify the malfunctioning part, four symbols are foreseen which are described in Chapter 9.1 on page 38.

EN

5.3 Safety Requirements Before and During Use

⚠ WARNING: Before starting to operate with the device, make sure that you have stock material (handpiece, inserts, wrenches) available to use in case the device stops due to a fault or of inconveniences.

⚠ CAUTION: Contraindication. Do not perform treatments on prosthetic artifacts made of metal or ceramics. The ultrasonic vibrations could lead to the de-cementing of the artifacts.

⚠ WARNING: Contraindications. Do not use the Piezosurgery® touch on patients who carry heart stimulators (Pace-makers) or other implantable electronic devices. This precaution also applies to the operator.

⚠ WARNING: Checking device status before the treatment. Always check that there is no water underneath the device. Before every treatment, always check that the device works perfectly and that the accessories are efficient. In case you uncover operating abnormalities, do not perform the treatment. Contact an Authorized Mectron Service Center if the abnormalities concern the device.

⚠ CAUTION: Clean function. After the device is used with aggressive and non-aggressive solutions, it is necessary to perform a cleaning cycle on the tubes and the handpiece with the clean function (see Cleaning and Sterilization Manual). If the tubes are not cleaned, the crystallization of the salts may seriously damage the device.

⚠ CAUTION: Contraindications. Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ WARNING: Infections control. First use: All new and repaired accessories are supplied in NON STERILE conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilization Manual. **Subsequent uses:** After every treatment, clean and sterilize all the reusable parts and accessories, following the instructions provided in the Cleaning and Sterilization Manual.

⚠ CAUTION: The electrical contacts inside the cord connector must be dry. Before connecting the handpiece to the device, make sure that the electrical contacts of the connector are perfectly dry, especially after the sterilization cycle in autoclave. If need be, dry the contacts by blowing compressed air onto them.

⚠ CAUTION: CLEAN function. The CLEAN function must be used after every treatment, before starting the cleaning and sterilization procedures.

⚠ **WARNING:** To ensure that the handpiece cools down, always activate it with the irrigation circuit correctly installed and filled. To fill the irrigation circuit, always use the PUMP function

⚠ **WARNING: Treatments that require irrigation.** Always check operation of the irrigation before and during use. Make sure the fluid outflows from the insert. Do not use the device if the irrigation does not work or if the pump is defective.

⚠ **WARNING: Only use original Mectron inserts, accessories, and spare parts.**

⚠ **WARNING: Use of non-original Mectron inserts:** this use entails finite damage to the handpiece threading, thus compromising correct operation and risking to cause harm to the patient.

⚠ **CAUTION:** To correctly use the device, you must press the foot pedal and start it with the insert not making contact with the part to be treated, so that the electronic circuit can recognize the best point of resonance of the insert without interferences, allowing its optimal performance.

⚠ **WARNING:** Before every treatment, make sure that the insert appropriate for the treatment is inserted on the handpiece. Exclusively use the Mectron torque wrench to fasten the insert to the handpiece. Do not use any other tools, such as pliers, pincers, etc.

⚠ **WARNING:** The patient must not come into contact with the device body or the foot pedal.

⚠ **WARNING:** Do not change the insert while the handpiece is operating, to prevent causing wounds to the operator.

⚠ **CAUTION:** Never force the cord-handpiece connector into its seat on the device as the cord-handpiece connector and/or the device can be damaged. If the connector cannot be easily fitted into its seat, probably the connector and its seat do not correctly match. Ensure that the connector index point of the cord-handpiece is placed upwards.

⚠ **CAUTION:** To avoid damaging the foot pedal cable, connect and/or disconnect it, holding it by the connector only. Never pull the cable.

⚠ **CAUTION:** Do not twist or rotate the foot pedal cable connector when inserting or removing it. Twisting can damage the connector.

⚠ **WARNING: Breakage and wear-out of the inserts.** High frequency oscillations and wear-out may, in rare circumstances, lead to the breakage of the insert. Do not bend, change shape of, or re-sharpen an insert in any way. Bending an insert or applying leverage on it can lead to its breakage. Deformed or otherwise damaged inserts are susceptible to breakage during their use. These inserts must never be used.

Excessive pressure on the inserts during their use can lead to their breakage. If an insert breaks, check that none of its fragments remain in the treated part and, at the same time, apply effective suction to remove them. The patient must be instructed to breathe through his nose during the treatment, or a dental dam must be used to prevent the patient from ingesting fragments of broken inserts. When the nitride coating wears out, the cutting efficiency decreases; re-sharpening the insert damages it and is therefore forbidden. Check that the insert is not worn out. Use of a worn-out insert reduces the cutting performance and can cause necrosis of the bone surface treated. During the intervention, frequently check that the insert is intact, especially its top. During the intervention, avoid prolonged contact with retractors or with metallic instrumentation in use.

⚠ **CAUTION: Contraindications.** Allow reusable, autoclavable items (the handpiece, the inserts, the torque wrench, and any other accessory that can be sterilized) to gradually return to room temperature after steam sterilization and prior to usage. The cooling process must not be accelerated.

⚠ WARNING: Check the level of physiological solution contained in the irrigation bag. Replace the irrigation bag with a new one before it is completely empty.

⚠ WARNING: Sterile single-use material. Before surgery. Before using any insert, check the sterile package integrity. Insert loses sterility in case of broken or damaged packaging.

⚠ WARNING: The PIEZOSURGERY® irrigation set is guaranteed for a single use, only. Separate and dispose of this kit according to the current regulations governing hospital waste.

⚠ WARNING: Before using the device, check that the irrigation tube clamps/roller are open. After surgery, close the clamps/roller before disconnecting the irrigation kit from the physiological irrigation bag.

⚠ WARNING: Be very careful to the sharp edges of cutting inserts. When tightening and removing these insert types, the edges may cause injuries.

⚠ CAUTION: Intermittent function. Prolonged use can cause the handpiece to overheat. Refer to *Chapter 8 on page 29* for average use times (intermittent function).

⚠ WARNING: Sterile single-use. Sterile single-use objects must only be used in a single surgical operation and for a patient, only. Single-use objects must not be reused. Separate and dispose of each single-use object according to current regulations governing hospital waste.

⚠ WARNING: During surgery, do not perform any maintenance activity on the system.

⚠ WARNING: The device is classified for intermittent mode operation (as reported in *Chapter 8 on page 29*). The continuous use of the device for a prolonged time and in any case exceeding the declared time limits can cause overheating, in particular, of the handpiece. In case of overheating, avoid contact of the handpiece with the operator and the patient.


5.4 Instructions for Use

After having connected all the accessories as described in Chapter 4.2 on page 12, proceed as follows:

EN


Open the air inlet on the spike, before proceeding with the operation.
Open the irrigation tube clamps/roller, if closed;

1


To fill the irrigation circuit, use the PUMP function by selecting PUMP/CLEAN on the touch keyboard: all the other selection options present on the display are disabled, and the PUMP/CLEAN writing flashes;

2


While the PUMP/CLEAN wording is flashing, press the foot pedal once and release it; the PUMP/CLEAN stops flashing and the irrigation circuit starts to fill up;

3


As soon as the peristaltic pump starts, the entire scale of values of the "irrigation" section lights up and during the liquid passage, the value of the irrigation shifts from 6 to 0;


4


5 The cycle can be interrupted as soon you see liquid out-flowing from the PIEZOSURGERY® touch handpiece, by pressing PUMP/CLEAN again or, alternatively, by pressing the foot pedal. The Pump function is disabled and the keyboard is enabled again, and displays the last setting used;


6 Screw the chosen insert onto the PIEZOSURGERY® touch handpiece, till the end of the threading;

7 Tighten the insert by using the Mectron torque wrench. To correctly use the Mectron torque wrench, operate as follows:

8 Place the insert into the wrench as shown;

Hold the central body of the handpiece firmly.


⚠ CAUTION: Do not grasp the end and/or cord of the handpiece; only hold the central body. Do not rotate the handpiece; simply hold it firmly and rotate the wrench, only.

Turn the wrench in a clockwise direction until the clutch engages (the outer part of the wrench slips, compared to the handpiece, producing mechanical a clicking sounds "CLICK"). The insert is now properly tightened;


On the keyboard, select the type of function, the irrigation level and the light option, if required.

⚠ CAUTION: For the correct operating parameters to be adopted for a specific insert, refer to the Chart annexed to this manual titled "Appropriate settings for the inserts on the PIEZOSURGERY® touch" or the illustrative leaflet of the purchased Mectron insert.

9

10

EN

5.5 Important Information on Inserts

WARNING:

- The insert must always be kept moving. If the insert is blocked or kept in contact for a long time, thermal damage of the treated part can be generated (for example, in the alveolar pockets, in the interdental septa or in the alveolar walls during the extraction procedures). We recommend using a continuous movement to minimise contact between the tip and the tissue. We recommend using high levels of irrigation as the power level is increased.
- For best performances, apply moderate, constant pressure on the insert. Do not apply excessive pressure; let the ultrasonic vibrations operate.
- When the titanium nitride coating is visibly worn out, the insert must be replaced. Use of an overly worn out insert reduces its cutting efficiency.
- Diamond-coated inserts: the diamond inserts must be replaced when the titanium nitride coating is visibly worn out and, at any rate, after a maximum of 10 treatments.
- Do not activate the handpiece while the insert is in contact with the treated part, so that the electronic circuit can recognize the best point of resonance of the insert and allow its optimal performance.
- Check the condition of wear of the insert and that it is intact before and during every use. In the event that a drop in performance becomes apparent, proceed to replace it.
- Use original Mectron inserts only. Use of non-original inserts, in addition to

voiding the warranty, damages the threading of the PIEZOSURGERY® touch handpiece, with the risk of no longer being able to screw the original inserts correctly during subsequent use. Moreover, the device settings are tested and guaranteed to operate correctly only when original Mectron inserts are used.

- Do not change the shape of the insert in any way by either bending or filing it. This could cause it to break.
- Do not use an insert that has suffered any type of deformation.
- Do not attempt to sharpen the insert used.
- Always check that the threaded parts of the insert and of the handpiece are perfectly clean – see the Cleaning and Sterilization Manual.
- If excessive pressure is applied on the insert, it can cause the insert to break and possibly harm the patient.
- For information on how to correctly use the inserts, refer to the annexed sheet “Appropriate settings for the inserts on the PIEZOSURGERY® touch” or the illustrative leaflet of the Mectron insert you’ve purchased.
- Before using the PIEZOSURGERY® touch, make sure you have prepared the operatory site by having first moved away the soft tissues, to avoid damaging them. It may happen that, while cutting the bone, accidental contact of certain parts of the insert with the soft tissues inflicts small traumas. To minimize this risk, use specific protective instruments.

6 MAINTENANCE

If the device is not used for a long time, observe the following recommendations:

1. Run a complete cleaning cycle on the irrigation circuit with the CLEAN function (see the Cleaning and Sterilization Manual);
2. Disconnect the device from the mains;
3. In case of a long period of inactivity, place the device in its original packaging and in a safe area;
4. Before using the device again, clean and

sterilise the handpiece and the wrench according to the instructions provided in the Cleaning and Sterilisation Manual;

5. Check that the inserts are not worn out, deformed, or broken, placing special attention to the integrity of their tip.

⚠ WARNING: Regularly check that the power supply cord is in good condition. If damaged, replace it with an original Mectron spare part.

7 DISPOSAL METHOD AND PRECAUTIONS

⚠ WARNING: Hospital waste.
Treat the following items as hospital waste:

- Inserts: when they are worn out or broken;
- Irrigation set: at the end of each intervention;
- Peristaltic pump tube: after 8 sterilization cycles;
- Torque wrench for inserts: when worn out or broken.

Disposable materials and materials that imply a biological risk  must be disposed of in

accordance with current local regulations governing hospital waste.

The PIEZOSURGERY® touch must be disposed of and treated as a waste requiring separate collection.

The customer is entitled to deliver the old device for disposal to the retailer supplying the new equipment. Instructions on correct disposal are available from Mectron.

Non-compliance with the previous points may produce a fine in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive.

8 TECHNICAL SPECIFICATIONS

| | |
|---|--|
| Device compliant with Regulation (EU) 2017/745 | Classe IIa |
| Classification as per IEC/EN 60601-1 | I Applied part type B (insert) IP 20 (device) IP 22 (foot pedal - FS-01 model) |
| Essential performance | According to the standard IEC 80601-2-60 the device has no essential performance. |
| Device for intermittent operation | 60sec. ON - 30sec. OFF with irrigation 30sec. ON - 120sec. OFF without irrigation (ENDO, PERIO) |
| Power supply voltage | 100-240 V~ 50/60 Hz |
| Max. power consumption | 120 VA |
| Fuses | Type 5 x 20 mm, T 2AL, 250V |
| Operating frequency | Automatic scan From 24 KHz to 36 KHz |
| Power settings | ENDO PERIO SPECIAL CANCELLOUS CORTICAL IMPLANT |
| Peristaltic pump capacity | Adjustable on the touch screen: ENDO / PERIO - 7 flow levels: from 0 to 6 (from 0 to approximately 75ml/min) SPECIAL / CANCELLOUS / CORTICAL / IMPLANT - 6 flow levels: from 1 to 6 (from 8 to approximately 75ml/min) |
| LED system of the handpiece | Light function on AUTO: LED system of the handpiece switches on as soon as the device starts to operate and switches off 3 seconds after the foot pedal is released. Light function on ON: The LED of the handpiece is always on; after 100 seconds of foot pedal disuse it switches off by itself and the light function moves into position AUTO. Light Function on OFF: The LED of the handpiece is always off. White LED light power risk free according to IEC/EN 62471. |
| APC circuit protections | No handpiece detected Cord interruption Insert not tightened correctly or broken |

| | |
|---|---|
| Operating conditions | from 10 °C to 40 °C Relative humidity from 30% to 75% Air pressure P: 800hPa/1060hPa |
| Transport and storage conditions | from -10 °C to 60 °C Relative humidity from 10% to 90% Air pressure P: 500hPa/1060hPa |
| Peristaltic pump tubing | It is advisable not to exceed 8 sterilization cycles |
| Altitude | Less than or equal to 2000 meters |
| Weight and size | 3,2Kg 300 x 250 x 95 mm (L x l x H) ^{a)} |

Table 2 – Technical Specifications

a) W = width ; L = length ; H = height

8.1 Electromagnetic Compatibility IEC/EN 60601-1-2

⚠️ WARNING: Contraindications. Interference with other equipment.

Although PIEZOSURGERY® touch complies with standard IEC 60601-1-2, it may interfere with other equipment located nearby. PIEZOSURGERY® touch must not be used near, or stacked on other equipment. Install PIEZOSURGERY® touch at a safe distance from life-support equipment. If necessary, however, prior their use, please check and monitor that the device and all other equipment function correctly in that layout.

⚠️ WARNING: Portable and mobile radio communication appliances may affect proper device operation.

⚠️ WARNING: Contraindications. Interference from other equipment.

An electrical scalpel or other electro-surgical units near the PIEZOSURGERY® touch device may interfere with its correct operation.

⚠️ WARNING: The device requires specific EMC precautions and must be installed and commissioned according to the EMC information provided in this chapter.

⚠️ WARNING: The use of cables and accessories not supplied by Mectron might negatively affect the EMC performances.

8.1.1 Guidance and Manufacturer’s Declaration - Electromagnetic Emissions

PIEZOSURGERY® touch is intended for use in the electromagnetic environment specified below. The customer or the user of PIEZOSURGERY® touch should ensure that it is used in such an environment.

| Emissions Test | Compliance | Electromagnetic Environment Guidance |
|---|------------|--|
| RF Emissions CISPR 11 | Group 1 | PIEZOSURGERY® touch only uses RF energy for internal function. Therefore, its FR emissions are very low and are not likely to cause any interference with nearby electronic equipment. |
| RF Emissions CISPR 11 | Class B | PIEZOSURGERY® touch is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Emissions of fluctuations voltage/flicker IEC 61000-3-3 | Compliant | |

8.1.2 Accessible Parts of the Casing

PIEZOSURGERY® touch is designed to operate in the electromagnetic environment specified below.

The customer or the user of PIEZOSURGERY® touch should ensure that it is used in such an environment.

| Phenomenon | Basic EMC standard or test method | Immunity test levels | Electromagnetic Environment Guidance |
|--|-----------------------------------|--|--|
| Electrostatic discharge (ESD) | IEC 61000-4-2 | ±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air | The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30 %. |
| Radiated RF EM fields ^{a)} | IEC 61000-4-3 | 3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80 % AM a 1 kHz ^{c)} | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |
| RATED power frequency magnetic fields ^{d) e)} | IEC 61000-4-8 | 30 A/m ^{g)} 50 Hz o 60 Hz | The magnetic fields at the mains frequency should have levels characteristic of a typical location in a commercial or hospital environment. |

a) The interface between the PATIENT physiological signal simulation, if used, and the PIEZOSURGERY® touch shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the PIEZOSURGERY® touch.

b) PIEZOSURGERY® touch that intentionally receives RF electromagnetic energy for the purpose of its operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.

c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

e) During the test, the PIEZOSURGERY® touch may be powered at any NOMINAL input voltage, but with the same frequency as the test signal.

f) Before modulation is applied.

g) This test level assumes a minimum distance between the PIEZOSURGERY® touch and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that PIEZOSURGERY® touch will be used closer than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.

8.1.3 Guidance and Manufacturer’s Declaration - Electromagnetic Immunity

8.1.3.1 Power Connection A.C. Input

PIEZOSURGERY® touch is designed to operate in the electromagnetic environment specified below.

The customer or the user of PIEZOSURGERY® touch should ensure that it is used in such an environment.

| Phenomenon | Basic EMC standard or test method | Immunity test levels | Electromagnetic Environment Guidance |
|---|-----------------------------------|--|--|
| Electrical fast transients / bursts ^{a) l) o)} | IEC 61000-4-4 | ±2 kV on contact 100 KHz repetition frequency | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Surges ^{a) b) j) o)} Line-to-line | IEC 61000-4-5 | ± 0.5 kV, ± 1 kV | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Surges ^{a) b) j) k) o)} Line-to-ground | IEC 61000-4-5 | ± 0.5 kV, ± 1 kV, ± 2kV | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Conducted disturbances induced by RF fields ^{c) d) o)} | IEC 61000-4-6 | 3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in the ISM bands between 0.15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 KHz ^{e)} | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |
| Voltage dips ^{f) p) r)} | IEC 61000-4-11 | 0% UT; 0,5 cycle ^{g)} At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle ^{e)} 70 % UT; 25/30 cycle ^{h)} Single phase: at 0° | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Voltage interruptions ^{f) i) o) r)} | IEC 61000-4-11 | 0% UT; 250/300 cycle ^{h)} | The quality of the network voltage should be that of a typical commercial or hospital environment. |

a) The test may be performed at any one power input voltage within the PIEZOSURGERY® touch RATED voltage range. If the PIEZOSURGERY® touch is tested at one power input voltage, it is not necessary to re-test at additional voltages.

b) All PIEZOSURGERY® touch cables are attached during the test.

c) Calibration for current injection clamps shall be performed in a 150 Ω system.

- d) If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) ME EQUIPMENT and ME SYSTEMS with a d.c. power input intended for use with a.c.-to-d.c. converters shall be tested using a converter that meets the specifications of the MANUFACTURER of the ME EQUIPMENT or ME SYSTEMS. The IMMUNITY TEST LEVELS are applied to the a.c. power input of the converter.
- g) Applicable only to PIEZOSURGERY® touch connected to single-phase a.c. mains.
- h) E.g. 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). PIEZOSURGERY® touch with battery backup shall resume line power operation after the test. For ME EQUIPMENT and ME SYSTEMS with RATED input current not exceeding 16 A, all phases shall be interrupted simultaneously.
- j) ME EQUIPMENT and ME SYSTEMS that does not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and ± 1 kV line(s) to line(s).
- k) Not applicable to CLASS II PIEZOSURGERY® touch.
- l) Direct coupling shall be used.
- m) r.m.s. , before modulation is applied.
- n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- o) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.
- p) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase.
- q) At some phase angles, applying this test to ME EQUIPMENT with transformer mains power input might cause an overcurrent protection device to open. This can occur due to magnetic flux saturation of the transformer core after the voltage dip. If this occurs, the ME EQUIPMENT shall provide BASIC SAFETY during and after the test.
- r) For ME EQUIPMENT and ME SYSTEMS that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the minimum and maximum RATED input voltage. ME EQUIPMENT and ME SYSTEMS with a RATED input voltage range of less than 25 % of the highest RATED input voltage shall be tested at one RATED input voltage within the range.

8.1.3.2 Points of Contact with the Patient

PIEZOSURGERY® touch is designed to operate in the electromagnetic environment specified below.

The customer or the user of PIEZOSURGERY® touch should ensure that it is used in such an environment.

| Phenomenon | Basic EMC standard or test method | Immunity test levels | Electromagnetic Environment Guidance |
|---|-----------------------------------|--|--|
| Electrostatic discharges (ESD) ^{c)} | IEC 61000-4-2 | ±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air | The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30 %. |
| Conducted disturbances induced by RF fields ^{a)} | IEC 61000-4-6 | 3 V ^{b)} 0.15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 KHz | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |

a) The following apply:

- All PATIENT-COUPLED cables shall be tested, either individually or bundled
- PATIENT-COUPLED cables shall be tested using a current clamp unless a current clamp is not suitable. In cases where a current clamp is not suitable, an EM clamp shall be used.
- No intentional decoupling device shall be used between the injection point and the PATIENT COUPLING POINT in any case.
- Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- Tubes that are intentionally filled with conductive liquids and intended to be connected to a PATIENT shall be considered to be PATIENT-COUPLED cables.

- If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18, 17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29, 7 MHz and 50,0 MHz to 54,0 MHz.

- b) R.M.S., before modulation is applied.
- c) Discharges shall be applied with no connection to an artificial hand and no connection to PATIENT simulation. PATIENT simulation may be connected after the test as needed in order to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

8.1.3.3 Parts Accessible to the Input / Output Signals

PIEZOSURGERY® touch is designed to operate in the electromagnetic environment specified below.

The customer or the user of PIEZOSURGERY® touch should ensure that it is used in such an environment.

| Phenomenon | Basic EMC standard or test method | Immunity test levels | Electromagnetic Environment Guidance |
|---|-----------------------------------|---|--|
| Electrostatic discharges (ESD) ^{e)} | IEC 61000-4-2 | ±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air | The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30%. |
| Electrical fast transients / bursts ^{b) f)} | IEC 61000-4-4 | ±1 kV on contact 100 KHz repetition frequency | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Surges Line-to-ground ^{a)} | IEC 61000-4-5 | ± 2kV | The quality of the network voltage should be that of a typical commercial or hospital environment. |
| Conducted disturbances induced by RF fields ^{b) d) g)} | IEC 61000-4-6 | 3 V ^{h)} 0.15 MHz - 80 MHz 6 V ^{h)} in the ISM bands between 0.15 MHz and 80 MHz ⁱ⁾ 80 % AM a 1 KHz ^{c)} | Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter. |

- a) This test applies only to output lines intended to connect directly to outdoor cables.
- b) SIP/SOPS whose maximum cable length is less than 3 m in length are excluded.
- c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) Calibration for current injection clamps shall be performed in a 150 Ω system.
- e) Connectors shall be tested per 8.3.2 and Table 4 of IEC 61000-4-2:2008. For insulated connector shells, perform air discharge testing to the connector shell and the pins using the rounded tip finger of the ESD generator, with the exception that the only connector pins that are tested are those that can be contacted or touched, under conditions of INTENDED USE, by the standard test finger shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) Capacitive coupling shall be used.

- g) If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- h) R.M.S., before modulation os applied.
- i) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18, 17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29, 7 MHz and 50,0 MHz to 54,0 MHz.

8.1.4 Specifications of the tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

PIEZOSURGERY® touch is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The customer or the user of PIEZOSURGERY® touch can help prevent electromagnetic interference by ensuring a minimum distance between the mobile and portable RF (transmitters) communication devices and PIEZOSURGERY® touch, as recommended, in relation to the maximum output power of radiocommunications equipment.

| Test Freq. (MHz) | Band ^{a)} (MHz) | Service ^{a)} | Modulation ^{b)} | Max power (W) | Distance (m) | Immunity test level (V/m) |
|------------------|--------------------------|--|---|---------------|--------------|---------------------------|
| 385 | 380 - 390 | TETRA 400 | Pulse modulation ^{b)} 18 Hz | 1.8 | 0.3 | 27 |
| 450 | 430 - 470 | GMRS 460 FRS 460 | FM ^{c)} ± 5 kHz deviation 1 kHz sine | 2 | 0.3 | 28 |
| 710 | 704 - 787 | LTE band 13, 17 | Pulse modulation ^{b)} 217 Hz | 0.2 | 0.3 | 9 |
| 745 | | | | | | |
| 780 | | | | | | |
| 810 | 800 - 960 | GSM 800/900 TETRA 800 iDEN 820 CDMA 850 Band LTE 5 | Pulse modulation ^{b)} 18 Hz | 2 | 0.3 | 28 |
| 870 | | | | | | |
| 930 | | | | | | |
| 1720 | 1700 - 1990 | GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS | Pulse modulation ^{b)} 217 Hz | 2 | 0.3 | 28 |
| 1845 | | | | | | |
| 1970 | | | | | | |
| 2450 | 2400 - 2750 | Bluetooth WLAN 802.11 b/g/n RFID 2450 Band LTE 7 | Pulse modulation ^{b)} 217 Hz | 2 | 0.3 | 28 |
| 5420 | 5100 - 5800 | WLAN 802.11 a/n | Pulse modulation ^{b)} 217 Hz | 0.2 | 0.3 | 9 |
| 5500 | | | | | | |
| 5785 | | | | | | |

a) For some services, only the uplink frequencies are included.

does not represent actual modulation, it would be worst case.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the PIEZOSURGERY® touch may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

⚠ WARNING: Portable RF communication equipment (including peripheral devices as antenna cables and external antennas) must not be used closer than 30 cm (12 inches) to any part of the device PIEZOSURGERY® touch, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these devices.

EN

9 TROUBLESHOOTING

9.1 Diagnostic System and Symbols on Keyboard

PIEZOSURGERY® touch is equipped with a diagnostic circuit that allows to detect operating abnormalities and to view their type on the keyboard via their relative symbol. By

using the following chart, the user is guided toward the identification and possible solution of the malfunction detected.

| Symbol on the Keyboard | Possible Cause | Solution |
|---|--|--|
|  | Handpiece electric contacts/ lead wet. | Thoroughly dry the contacts with compressed air. |
| | PIEZOSURGERY® handpiece not connected to the device. | Connect the handpiece. |
| | Faulty Handpiece. | Replace the handpiece. |
| | Malfunction of the tuning circuit. | Contact an Authorized MECTRON Service Centre. |
|  | Insert not properly tightened on the handpiece. | Unscrew the insert and screw it back correctly using the torque wrench (see <i>Chapter 5.4 on page 24</i>). |
| | Insert broken, worn or deformed. | Replace the insert. |
| | Handpiece electric contacts/ lead wet. | Thoroughly dry the contacts with compressed air. |

| Symbol on the Keyboard | Possible Cause | Solution |
|---|--|---|
|  | Peristaltic pump malfunction. | Check that there are no impediments to pump rotation. |
| | Tubing not correctly positioned in the peristaltic pump. | Correctly reposition the silicon tubing inside the pump (see <i>Chapter 4.2</i> on page 12). |
| | The device has been turned off and on again without waiting 5 seconds. | Turn device off and wait 5 seconds before turning it on again. |
| | Abnormalities in the power distribution network, excessive electrostatic discharges or internal anomalies. | Turn the device off and wait 5 seconds before turning it on again. If the message persists, contact an Authorized Mectron Service Center. |
|  | The system was started with the pedal pressed. | Check that during system start-up the pedal is not pressed. If the problem persists, disconnect the pedal and contact an Authorized Mectron Service Center. |

Tabella 3 – Diagnostic System and Symbols on Keyboard

NOTE: For any error message not listed in the table, contact a Mectron Service Centre.

9.2 Troubleshooting Quick Guide

EN

| Problem | Possible Cause | Solution |
|--|--|--|
| The device does not turn on after having brought the switch into position "I" | The power supply connector is not properly plugged into the socket on the back of the device. | Check that the power supply cord is firmly connected. |
| | The power supply cord is faulty. | Check that the power supply socket is functioning. Replace the power supply cord. |
| | The fuses are out-of-order. | Replace the fuses (refer to <i>Chapter 9.3 on page 42</i>) |
| The device is switched on but not working. The display does not show any error. | The foot pedal plug is not correctly inserted into the device socket. | Correctly insert the foot pedal plug into the socket on the back of the device. |
| | The foot pedal is not working properly. | Contact the Authorised Mectron Service Centre. |
| The device is switched on but not working. One of the following symbols appears on the screen:  | See <i>Chapter 9.1 on page 38</i> for the possible cause, according to the symbol that has been displayed. | See <i>Chapter 9.1 on page 38</i> for the action to undertake, according to the symbol that has been displayed |
| When operating, the PIEZOSURGERY® handpiece emits a slight whistle. | The insert is not properly tightened into the handpiece. | Unscrew the insert and screw it back correctly by using the Mectron torque wrench. (Refer to <i>Chapter 5.4 on page 24</i>) |
| | The irrigation circuit has not been completely filled. | Fill the irrigation circuit by using the "pump" function (Refer to <i>Chapter 5.4 on page 24</i>) |

| Problem | Possible Cause | Solution |
|---|--|---|
| During operation, no liquid flows out from the insert. | The insert is of the type which does not provide the passage of liquids. | Use an insert of the type with the passage of liquids. |
| | The insert is clogged. | Unscrew the insert from the handpiece and unclog the insert water passage by blowing compressed air through it. If the problem persists, replace the insert with a new one. |
| | The handpiece is clogged. | Contact an Authorised Mectron Service Centre. |
| | The irrigation level on the screen is set to "0". | Adjust the irrigation level |
| | The liquid irrigation bag is empty | Replace the empty bag with a new one full. |
| | The air inlet on the spike/roller has not been opened. | Open the air inlet on the spike/roller. |
| | The silicon tubing are incorrectly installed. | Check the irrigation kit connections. |
| The device is working properly but the pump works with difficulty. | Excessive pressure of the peristaltic pump impeller on the tube. | Check that the tube in the peristaltic pump has been correctly inserted. (Refer to <i>Chapter 4.2 on page 12</i>) |
| The pump is working properly, but when it stops, liquid flows out from the handpiece. | The peristaltic pump cover is not closed properly. | Check that the peristaltic pump cover is closed properly. (Refer to <i>Chapter 4.2 on page 12</i>) |
| Poor performance | The insert is not correctly tightened on the handpiece | Unscrew the insert and screw it back correctly by using the Mectron torque wrench. (Refer to <i>Chapter 5.4 on page 24</i>) |
| | Broken, worn-out or deformed insert. | Replace the insert with a new one. |

Tabella 4 – Troubleshooting Quick Guide


9.3 Fuses Replacement

⚠ WARNING: Switch the device off.
Always turn the device off by means of the main switch and disconnect it from the electrical power socket before proceeding.

EN

Apply leverage with a flat screwdriver, inserting its tip in the seat of the fuse-holder drawer located under the power supply socket;


1

Pull out the fuse-holder drawer;


⚠ WARNING: Replace the fuses with others respecting the characteristics indicated in Chapter 8 on page 29.

2

Reinsert the drawer in its housing.

3

9.4 Sending to an Authorized Mectron Service Centre

In the event that technical assistance is needed on the device, contact one of the Authorized Mectron Service Centers or your Retailer. Do not attempt to repair or modify the device and its accessories.

Clean and sterilize all the parts that need to be sent to an Authorized Mectron Service Center, following the instructions provided in the Cleaning and Sterilization Manual supplied with the device.

Leave the sterilized parts in the packet, procedure that confirms the sterilization process has been performed.

The cleaning and sterilization requirements are in line with those in force concerning the safeguarding of health and safety on the workplace, as per DLgs 81/08 and its subsequent amendments, both laws of the Italian State.

In the case that the customer does not fulfill the requirements, Mectron reserves the right to charge him or her the cleaning and sterilization expenses, or to reject the goods received in unsuitable conditions, returning them to the customer, at his or her expense, for them to be correctly cleaned and sterilized. The device must be returned properly packaged, accompanied by all the accessories

and by a sheet bearing the:

- Data of the owner with telephone number;
- Product name;
- Serial number and/or lot number;
- Reason for goods returned / description of the malfunction;
- Photocopy of delivery note or purchase invoice of the device.

⚠ CAUTION: Packaging

Pack the device in its original package to prevent damages during transport.

Once the material has been received by the Authorized Mectron Service Center, qualified technical personnel will provide its evaluation to the given circumstance. The repair works will be performed only upon their prior acceptance by the end customer. For further details, contact the closest Authorized Mectron Service Center or your retailer.

Unauthorized repair works may damage the system and make the warranty void, and absolve Mectron of any responsibility for direct or indirect harm to persons or damages to things.

10 WARRANTY

All the Mectron devices are subjected to an accurate end inspection that ascertains their full functionality before they are placed on the market.

Mectron provides warranty on PIEZOSURGERY® touch purchased new from a retailer or a Mectron importer to cover material or manufacturing defects for:

- 2 YEARS (TWO) on the device from the date of purchase;
- 1 YEAR (ONE) on the handpiece complete with cord from the date of purchase.

The accessories are not included in the warranty.

During the period in which the warranty is valid, Mectron commits itself to repair (or at its own free choice, replace) those product parts that were to prove defective according to its judgment, free of charge.

Entire replacement of the Mectron products is excluded.

The warranty of the manufacturer and the homologation of the device are not valid in the following cases:

- The device is not used in accordance to the intended use foreseen for it.
- The device is not used in compliance to all the instructions and precautions described in this manual.
- The electrical system in the premises in which the device is used is not compliant to the norms in force and to the relative precautions.
- The assembly operations, extensions, adjustments, updates, and repair works are performed by personnel not authorized by Mectron.

- The environmental conditions for the preservation and storage of the device are not compliant to the precautions indicated in the Chapter 8 on page 29.
- Use of non-original Mectron inserts, accessories, and spare parts that can compromise correct device operation and cause harm to the patient.
- Accidental breakage due to transport.
- Damages due to incorrect use or negligence, or to connection to a power voltage other than the one foreseen.
- Warranty expiry.

The expected service life of the device is 5 years, minimum.

The service life/duration does not define a limit of use; the service life of the device defines the period of time, subsequent to installation and/or commissioning, during which the original performances or, in any case, adequate for the intended use are guaranteed, without any degradation occurring such as to compromise its functionality and reliability.

The service life is a minimum qualitative target of the design, therefore, it is not excluded that single parts or components guarantee performances and reliability higher than those declared by the manufacturer.

The service life is intended in compliance with the maintenance plans provided for in this manual, it does not include components normally subjected to "wear" and it is independent of the warranty period: the service life period does not establish any implicit or explicit extension of the warranty period.

CAUTION

The warranty starts from the date of purchase of the device, which evidence is given by the delivery note/purchase invoice issued by the Dealer / Importer or, in case of device with activation code, from the date of activation of the same.

In order to avail of the warranty service, the customer must return, at its own expense, the device to be repaired to the MECTRON Dealer / Importer from which they purchased the product. The device must be returned together with the original packaging, accompanied by all the accessories and by a form containing:

- The data of the owner and telephone number;
- The data of the Dealer / Importer;
- Photocopy of the delivery note/purchase invoice of the device by the owner where are reported the date, the name of the device and the serial number;
- Description of the failure.

The transport and the damage caused by transport are not covered by the warranty.

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© Mectron S.p.A. 2021. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.