Ultra Rapid Readout
Fluorescence System

Hole with special paper /
Orificio con papel especial
Plastic cap /
Tapa plástica
Plastic tube /
Tubo plástico
Identification label /
Etiqueta identificadora
Process indicator /
Indicador de proceso
Glass ampoule and growth medium /
Ampolla de vidrio y medio de cultivo
Filter /
Filtro
Spore carrier /
Portador de esporas

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Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection.

Terragene® and Bionova® are registered trademarks of Terragene S.A.

Lic. Adrián J. Povetto
Director Técnico
Technical Director

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Explanation of Symbols

Product is designed for use with Steam sterilization cycles.
 Batch number
 Manufacture Date.
 Expiration Date.
 Manufacturer.

Intended Use Table

Model	Autoclave/Steam Cycles	
	Gravity Displacement	
	132 °C	135 °C
	10 minutes	
	15 minutes	10 minutes
	25 minutes	
BT224	Dynamic Air Removal (Vacuum Assist)	
	132 °C	135 °C
	4 minutes	3 minutes
Fluorescence Read Time	pH Color Change	
	20 minutes	48 hours

EN Biological Indicators
For Steam sterilization

Indications for Use

United States

Terragene® Bionova® SCBI (BT224) is a self-contained biological indicator inoculated with viable 10⁶ *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of steam sterilization processes. On each Terragene® Bionova® SCBI is a chemical process indicator that changes color from pink to brown when exposed to Steam. See Intended Use Table for reference.

Outside the United States

Terragene® Bionova® BT224 Fluorescence Ultra Rapid Readout Biological Indicators have been designed for quick and easy monitoring of vacuum assisted and gravity displacement steam sterilization cycles at 132-135 °C.

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7. Colocar el Process Indicator on SCBI label. A color change to brown indicates that the SCBI has been exposed to Steam. **IMPORTANT:** This color change does not evidence the process effectiveness to achieve sterility. If the Process Indicator color has not changed, check the sterilization process.

8. Press the lid to seal the tube. Crush the ampoule contained in the SCBI with an individual ampoule crusher or the ampoule crusher placed within the incubator's incubation area. Then shake the tube down vigorously with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the SCBI in the incubator. **IMPORTANT:** Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that correct incubation conditions were met; capability of medium to promote rapid growth; viability of spores after being heated to 132 °C; and compatibility of the medium to chemicals and proper functioning of Terragene® Bionova® Auto-Reader Incubators. Both, the positive control indicator and the processed indicator, should belong to the same batch.

9. Incubar el indicador procesado y el indicador usado como control positivo en el appropriate Terragene® Bionova® Auto-Reader Incubators for a maximum of 20 minutes at (60 ± 2) °C for ultra rapid readout. **NOTE:** Holding time between sterilization and incubation should not exceed a 7-day period. Fluorescence detection by the reader (excitation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at 20-minute incubation, the result is negative. The positive control must give positive fluorescence readout. It is good practice to incubate a positive control for a visual color change.

10. Record the positive results and discard the SCBI immediately, as it is indicated below.

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Lectura Ultra Rápida: 20 minutos

The Ultra Rapid Readout must be carried out in the appropriate Terragene® Bionova® Auto-Reader Incubators. Fluorescence is emitted after the reader stimulates the spores with UV light. Final readout of negative results is readily available after 20 minutes of incubation. Fluorescence readout is an indirect measure of the germination and growth of *Geobacillus stearothermophilus* spores, which have resisted the sterilization process (positive result). Furthermore, a failure in the sterilization process can also become evident by culture medium color change. Due to the high sensitivity of the fluorescence results at 20 minutes, conventional incubation for color change is not an advantage.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn yellow during incubation at (60 ± 2) °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after incubation. The positive control must show a color change from purple to yellow for results to be valid.

Readout time: 7 days

A 7-day readout is optional and may be routinely performed. This is an initial validation of the 20-minute readout. Fluorescence results are compared to the 7-day visual readouts. **NOTE:** If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

Monitoring frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As the best practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate biological indicator.

Storage

Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI.

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